

PCR screening of MiMIC RMCE events to determine orientation of GFSTF or Trojan-GAL4 cassette

This is the protocol that we use to screen and verify the stocks we send to the Bloomington Drosophila Stock Center. For details of the insertion collection see:

<http://flypush.imgen.bcm.tmc.edu/pscreen/rmce/>

Venken et al. (2011) *Nature Methods* 8:737-743

Nagarkar-Jaiswal et al. (2015) *eLife* 4:e05338

Nagarkar-Jaiswal et al. (2015) *eLife* 4:e08469

We recommend screening at least 5 independent candidate events.

Step 1. Genomic DNA extraction

We use the PureLink® Genomic DNA Mini Kit from Thermo Fisher (cat. no. K1820-02 for 250 reactions). All solutions mentioned below (with the exception of EtOH) are provided in the kit.

1. Prepare proteinase K digestion reaction master mix (total 200 µl per reaction.)
 - a. Digestion Buffer - 180 µl per reaction (180 x ___ candidates)
 - b. Proteinase K - 20 µl per reaction (20 x ___ candidates)
2. Collect 10 to 15 adult flies for each candidate and homogenize with 200 µl digestion mix. We use polypropylene Kimble™ Kontes™ Pellet Pestle™ tissue grinders (Fisher Scientific cat. no. K749521-1590).
3. Incubate at 55°C for 60 min.
4. To each tube, add 20 µl RNase A. Stand at RT for 2 min.
5. Centrifuge at max for 5 min and transfer supernatant to fresh microcentrifuge tube.
6. To each tube, add:
 - a. Lysis/Binding Buffer 200 µl
 - b. 100% EtOH 200 µl
7. Vortex.
8. Run lysate through column (using centrifuge or vacuum manifold).
9. Wash with 500 µl Wash Buffer 1. (Spin at max for 1 min. Discard flow through.)
10. Wash with 500 µl Wash Buffer 2. (Spin at max for 1 min. Discard flow through.)
11. Spin once again at max for 1 min and discard flow through.
12. Place column in a collection tube, add 100 µl Elution Buffer and stand at RT for 1 min. Spin at max for 1 min.

Step 2. PCR reaction

1. Prepare PCR master mixes. (20 μ l reaction volume each candidate). We are using GoTaq® G2 Green Master Mix (Promega M7823, 1,000 reactions) for high throughput, but you can modify to use your lab's PCR Taq/dNTPs.

Reaction 1		
Template	1.0 μ l	Mix:
Primer 1	0.5 μ l	μ l
Primer 3 or 5	0.5 μ l	μ l
GoTaq Green	10.0 μ l	μ l
ddH ₂ O	8.0 μ l	μ l

Reaction 2		
Template	1.0 μ l	Mix:
Primer 1	0.5 μ l	μ l
Primer 4 or 6	0.5 μ l	μ l
GoTaq Green	10.0 μ l	μ l
ddH ₂ O	8.0 μ l	μ l

Reaction 3		
Template	1.0 μ l	Mix:
Primer 2	0.5 μ l	μ l
Primer 3 or 5	0.5 μ l	μ l
GoTaq Green	10.0 μ l	μ l
ddH ₂ O	8.0 μ l	μ l

Reaction 4		
Template	1.0 μ l	Mix:
Primer 2	0.5 μ l	μ l
Primer 4 or 6	0.5 μ l	μ l
GoTaq Green	10.0 μ l	μ l
ddH ₂ O	8.0 μ l	μ l

Primer sets:

Primer name	Sequence	Use for	Description
1 MiLF	GCGTAAGCTACCTTAATCTCAAGAAGAG	All	Forward primer for Minos
2 MiLR	CGCGGCGTAATGTGATTTACTATCATAC	All	Reverse primer for Minos
3 TagR	GTGGCTGTTGAAGTTGTACTC	GFSTF	Reverse primer for EGFP
4 TagF	GGATGACGGCACCTACAAGAC	GFSTF	Forward primer for EGFP
5 T2A-R	GGGTTCTCCTCCACATCGC	TG4	Reverse primer for TG4
6 GAL4-5F	AACTGTGCATCGTGACCATC	TG4	Forward primer for TG4

2. Cycling conditions are as follows (may need to be adjusted for your own Taq).

GFSTF insertions:

94°C	10 min	
94°C	30 sec	35 cycles
60°C	30 sec	
72°C	1 min	
72°C	10 min	
12°C	hold	

Trojan-GAL4 (TG4) insertions:

94°C	10 min	
94°C	30 sec	
56°C	30 sec	35 cycles
72°C	1 min	
72°C	10 min	
12°C	hold	

Step 3. Run PCR products on agarose gel

- If gene and original MiMIC insertion are in the same orientation (i.e., either both + or both -), then the successful RMCE should screen positive for PCR reactions 1 and 4 only.
- If gene and original MiMIC insertion are in opposite orientations, then the successful RMCE should screen positive for PCR reactions 2 and 3 only.

See Supplementary Figure 2, Venken et al. (2011) *Nature Methods* 8:737 for an example.